



Recombinant Human BET1 homolog (BET1)

Product Code	CSB-CF002670HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O15155
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MRRAGLGEGVPPGNYGNYGYANSYGYSACEEENERLTESLRSKVTAIKSLSIEI GHEVKTQNKLLAEMDSQFDSTTGFLGKTMGKLIKLSRGSQTKLLCYMMLFSLF VFFIIYWIIKLR
Source	in vitro E.coli expression system
Target Names	BET1
Protein Names	Recommended name: BET1 homolog Short name= hBET1 Alternative name(s): Golgi vesicular membrane-trafficking protein p18
Expression Region	1-118
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.