



Recombinant Human Receptor activity-modifying protein 1 (RAMP1)

Product Code	CSB-CF019304HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60894
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	CQEANYGALLRELCLTQFQVDMEAVGETLWCDWGRTIRSYRELADCTWHMA EKLGCFWPNAEVD RFFLAVHGRYFRSCPISGRAVRDPPGSILYPFIVVPITVTL VTALVWQSKRTEGIV
Source	in vitro E.coli expression system
Target Names	RAMP1
Protein Names	Recommended name: Receptor activity-modifying protein 1 Alternative name(s): Calcitonin-receptor-like receptor activity-modifying protein 1 Short name= CRLR activity-modifying protein 1
Expression Region	27-148
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.