



Recombinant Human Glycophorin-C (GYPC)

Product Code	CSB-CF010076HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04921
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MWSTRSPNSTAWPLSLEPDPGMASASTTMHTTTIAEPDPGMSGWPDGRMET STPTIMDIVVIAGVIAAVAIVLVSLLFVMLRYMYRHKGTYHTNEAKGTEFAESAD AALQGDPALQDAGDSSRKEYFI
Source	in vitro E.coli expression system
Target Names	GYPC
Protein Names	Recommended name: Glycophorin-C Alternative name(s): Glycoconnectin Glycophorin-D Short name= GPD Glycoprotein beta PAS-2' Sialoglycoprotein D CD_antigen= CD236
Expression Region	1-128
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.