



Recombinant Human Receptor activity-modifying protein 2 (RAMP2)

Product Code	CSB-CF019305HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60895
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	QPLPTTGTPGSEGGTVKKNYETAVQFCWNHYKDQMDPIEKDWCDWAMISRPY STLRDCLEHFAELFDLGFNPLAERIIFETHQIHFANCSLVQPTFSDPPEDVLLA MIIAPICLIPFLITLVVWRSKDSEAQA
Source	in vitro E.coli expression system
Target Names	RAMP2
Protein Names	Recommended name: Receptor activity-modifying protein 2 Alternative name(s): Calcitonin-receptor-like receptor activity-modifying protein 2 Short name= CRLR activity-modifying protein 2
Expression Region	43-175
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.