



Recombinant Human Probetacellulin (BTC)

Product Code	CSB-CF002853HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35070
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	DGNSTRSPETNGLLCGDPEENCAATTTQSKRKGHF SRCPKQYKHYCIKGRCR FVVAEQTPSCVCD EGYIGARCERVDLFYLRGDRGQILVICLIAMVVFIILVIGVC TCCHPLRKRRKRKKKEEEMETLGKDITPINEDIEETNIA
Source	in vitro E.coli expression system
Target Names	BTC
Protein Names	Recommended name: Probetacellulin Cleaved into the following chain: 1. Betacellulin Short name= 2. BTC
Expression Region	32-178
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	10xHis-SUMO-tag
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.