



Recombinant Rat Mitochondrial fission 1 protein (Fis1)

Product Code	CSB-CF008684RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P84817
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MEAVLNELVSVEDLKNFERKFQSEQAAGSVSKSTQFEYAWCLVRSKYND DIR RGIVLLELLPKGSKEEQRDYVFYLAVGNRYRLKEYEKALKYVVRGLLQTEPQNN QAKELERLIDKAMKKDGLVGMALVGGMALGVAGLAGLIGLAVSKSKS
Source	in vitro E.coli expression system
Target Names	Fis1
Protein Names	Recommended name: Mitochondrial fission 1 protein Alternative name(s): FIS1 homolog Short name= rFis1 Tetratricopeptide repeat protein 11 Short name= TPR repeat protein 11
Expression Region	1-152
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.