



Recombinant Rat T-cell surface glycoprotein CD3 delta chain (Cd3d)

Product Code	CSB-CF004930RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P19377
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	FKIEVVEYEDKVFVNCNTSIRHLDGSVERWLTKNKSLILGKGILDPRGMYMCNG TEELAKEVSTVQVYYRMCQNCVELDSATLAGVIITDLIATLLLALGVYCFAGHET GRLSGAVDTQVLLKNEQLYQPLRDRDDAQYSRLGGNWPRNKRS
Source	in vitro E.coli expression system
Target Names	Cd3d
Protein Names	Recommended name: T-cell surface glycoprotein CD3 delta chain Alternative name(s): T-cell receptor T3 delta chain CD_antigen= CD3d
Expression Region	22-173
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.