



Recombinant Bovine Transmembrane protein 190 (TMEM190)

Product Code	CSB-CF023780BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	E1BJD3
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	NGIQGFFYPWSCEGDVWDRESCGGQAAIENPNLCLRLRCCYRDGVCYHQRP DENMRRKHMWALGWTCGGLLFLITSICLFWWARRHDMRLRPWFLKKGCDLS RTVSLLSKDRTPSEKKTSPVGSIPPAAPTEGALDVSGGTEGEGTEGGEETEGG DEDD
Source	in vitro E.coli expression system
Target Names	TMEM190
Protein Names	Recommended name: Transmembrane protein 190
Expression Region	22-180
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.