



Recombinant Mouse T-cell surface glycoprotein CD3 gamma chain (Cd3g)

Product Code	CSB-CF004933MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11942
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	QTNKAKNLVQVDGSRGDGSVLLTCGLTDKTIKWLKDGSIISPLNATKNTWNLG NNAKDPRGTYQCQGAKETSNPLQVYYRMCENCIELNIGTISGFIFAEVISIFFLA LGVYLIAGQDGVQRASDKQTLQNEQLYQPLKDREYDQYSHLQGNQLRKK
Source	in vitro E.coli expression system
Target Names	Cd3g
Protein Names	Recommended name: T-cell surface glycoprotein CD3 gamma chain Alternative name(s): T-cell receptor T3 gamma chain CD_antigen= CD3g
Expression Region	23-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.