



Recombinant Rat Group XVI phospholipase A1/A2 (Pla2g16)

Product Code	CSB-CF018089RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53817
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MPIPEPKPGDLIEIFRPMYSHWAIYVGDGYVIHLAPPSEIPGAGAASIMSALTDK AIVKKELLRDVAGKDKYQVNNKHDKKEYTPLPLNKIIQRAEELVGQEVLYRLTSE NCEHFVNELRYGVPRSDQVRDAVKVATVTGVGLAALGLIGVMLSRNKKQKQ
Source	in vitro E.coli expression system
Target Names	Pla2g16
Protein Names	Recommended name: Group XVI phospholipase A1/A2 EC= 3.1.1.32 EC= 3.1.1.4 Alternative name(s): H-rev 107 protein HRAS-like suppressor 3
Expression Region	1-160
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.