



Recombinant Mouse Platelet glycoprotein IX (Gp9)

Product Code	CSB-CF009690MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88186
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	TQACPRPCTCQSLETMGLKVNCEGQGLTALPVIPAHTRQLLLANNSLRVPPG AFDHLPLQLWDLVDVTHNPWHCDCSLTYLRLWLEDHMPEALMHVYCASPDLAT RRPLGQLTGYELGSCGWKLPPSWAYPGVWWDVSLVAVAVLGLLLAGLLNTF TESRN
Source	in vitro E.coli expression system
Target Names	Gp9
Protein Names	Recommended name: Platelet glycoprotein IX Short name= GP-IX Short name= GPIX Alternative name(s): Glycoprotein 9 CD_antigen= CD42a
Expression Region	17-177
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.