



Recombinant Sheep T-cell surface glycoprotein CD3 epsilon chain (CD3E)

Product Code	CSB-CF004931SH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29328
Product Type	Transmembrane Protein
Immunogen Species	Ovis aries (Sheep)
Sequence	QDDTEQNPYEVSISGNSVELTCPKDFENGIQWKRNNEQMKGHNKEYLLLDQF SEMESSGYQCLATEGNTEAAHTLYLKARVCKNCMEVNLLEVATIIIVVDICVTL GLLLLVYYWSKSRKAKATPMTRGAGAGGRPRGQNRERPPPVPNPDYEPIRK GQRDLYSGLNQRGV
Source	in vitro E.coli expression system
Target Names	CD3E
Protein Names	Recommended name: T-cell surface glycoprotein CD3 epsilon chain Alternative name(s): CD_antigen= CD3e
Expression Region	22-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.