



Recombinant Human Signal peptidase complex catalytic subunit SEC11A (SEC11A)

Product Code	CSB-CF020931HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P67812
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLSLDFLDDVRRMNKRQLYYQVLNFGMIVSSALMIWKGLMVTGSESPIVVLS GSMEPAFHRGDLLFLTNRVEDPIRVGEIVVFRIEGREIPIVHRVLKIHEKQNGHIK FLTKGDNNVDDRGLYKQGQHWLEKKDVVGRARGFVPYIGIVTILMNDYPKFK YAVLFLGLFVLVHRE
Source	in vitro E.coli expression system
Target Names	SEC11A
Protein Names	Recommended name: Signal peptidase complex catalytic subunit SEC11A EC= 3.4.21.89 Alternative name(s): Endopeptidase SP18 Microsomal signal peptidase 18 kDa subunit Short name= SPase 18 kDa subunit SEC11 homolog A SEC11-like protein 1 SP
Expression Region	1-179
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.