



Recombinant Mouse Natural killer cells antigen CD94 (Klrd1)

Product Code	CSB-CF012470MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O54707
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAVSRITRWRLMSVIFGIKCLFLMVTLGVLLINSFTIQNIQSTPSPTTTVEFQEVS ECCVCLDKWVGHCNCYFISKEEKSWKRSRDFCASQNSSLLQPQSRNELSF MNFSQTFWIGMHYSEKRNAWLWEDGTVPSKDLFPEFSVIRPEHCIVYSPSKS VSAESCENKNRYICKKLPI
Source	in vitro E.coli expression system
Target Names	Klrd1
Protein Names	Recommended name: Natural killer cells antigen CD94 Alternative name(s): Killer cell lectin-like receptor subfamily D member 1 CD_antigen= CD94
Expression Region	1-179
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.