



Recombinant Mouse Platelet glycoprotein Ib beta chain (Gp1bb)

Product Code	CSB-CF009686MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56400
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	PAPCSCAGTLVDCGRRGLTWASLPAAFPPDTTELVLGTGNNLTALPPGLLDALP ALRAAHLGANPWRCDRCRLPLRAWLAGRPERAPYRDLRCVAPPALRGRLLPY VAEDELRAACAPGLLCWGALVAQLALLVLGLLHALLLALLLGRRLRRLRARARAR SIQEFSLTAPLVAESARGGAS
Source	in vitro E.coli expression system
Target Names	Gp1bb
Protein Names	Recommended name: Platelet glycoprotein Ib beta chain Short name= GP-Ib beta Short name= GPIb-beta Short name= GPIbB Alternative name(s): CD_antigen= CD42c
Expression Region	27-206
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.