



Recombinant Mouse Pre T-cell antigen receptor alpha (Ptcra)

Product Code	CSB-CF018961MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C6B2
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	LPSGIAGTPFP SLAPPITLLVDGRQHMLVVCLVLDAAPPGLDNPVWFSAGNGS ALDAFTYGP SLAPDGTWTS LAQLSLPSEELEAWEPLVCHTRPGAGGQNRSTH PLQLSGESSTARSCFPEPLGGTQRQVLWLSLLRLLLFKLLLLLDVLLTCSHLRLH VLAGQHLQPPPSRKSLPPTHRIWT
Source	in vitro E.coli expression system
Target Names	Ptcra
Protein Names	Recommended name: Pre T-cell antigen receptor alpha Short name= pT-alpha Short name= pTa Alternative name(s): gp33 pT-alpha-TCR
Expression Region	17-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.