



Recombinant Rat Syndecan-2 (Sdc2)

Product Code	CSB-CF020889RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34900
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	ETRAELTSDKMYLDSSSIEEASGLYPIDDDDYSSASGSGAYEDKGSPDLTTS QLIPRISLTSAAPEVETMTLKTQSITPTQTESPEETDKKEFEISEAEEKQDPAVK STDVYTEKHSDNLFKRTEVLAAVIAGGVIGFLFAIFLILLLLIVYRMRKKDEGSYDL GERKPSSAAYQKAPTKEFYA
Source	in vitro E.coli expression system
Target Names	Sdc2
Protein Names	Recommended name: Syndecan-2 Short name= SYND2 Alternative name(s): Fibroglycan Heparan sulfate proteoglycan core protein Short name= HSPG CD_antigen= CD362
Expression Region	19-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.