



Recombinant Human Tumor necrosis factor ligand superfamily member 4 (TNFSF4)

Product Code	CSB-CF023994HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P23510
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MERVQPLEENVGNAARPRFERNKLLLVASVIQGLGLLLCFTYICLHFSALQVSH RYPRIQSIKVQFTEYKKEKGFILTSQKEDEIMKVQNNSVIINCDFYLISLKGYS QEVNISLHYQKDEEPLFQLKKVRSVNSLMVASLTYKDKVYLVNVTDDNTSLDDF HVNGGELILIHQNPGEFCVL
Source	in vitro E.coli expression system
Target Names	TNFSF4
Protein Names	Recommended name: Tumor necrosis factor ligand superfamily member 4 Alternative name(s): Glycoprotein Gp34 OX40 ligand Short name= OX40L TAX transcriptionally-activated glycoprotein 1 CD_antigen= CD252
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.