



Recombinant Rat Proheparin-binding EGF-like growth factor (Hbegf)

Product Code	CSB-CF010154RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06175
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	ESLERLRRGLAAATSNPDPPTGTTNQLLPTGADRAQEVQDLEGTDLDFKVAF SSKPQALATPGKEKNGKKRKGKGLGKKRDPCLKKYKDYCIHGECRYLKELRI PSCHCLPGYHGQRCHGLTLPVENPLYTYDHTTVLAVVAVVLSSVCLLVIVGLL MFRYHRRGGYDLESEEKVKLGMASSH
Source	in vitro E.coli expression system
Target Names	Hbegf
Protein Names	Recommended name: Proheparin-binding EGF-like growth factor Cleaved into the following chain: 1. Heparin-binding EGF-like growth factor Short name= 2. HB-EGF Short name= 3. HBEGF
Expression Region	24-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.