



Recombinant Mouse T-cell antigen CD7 (Cd7)

Product Code	CSB-CF004953MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P50283
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	QDVHQSPRLTIASEGDSVNITCSTRGHLEGILMKKIWPQAYNVIYFEDRQEPTV DRTFSGRINFSGSQKNLTITISSLQLADTGDYTCFAVRKVSARGLFTTVVVKKEK SSQEAYRSQEPLQTSFSPAAIAVGFFFTGLLLGVVCSMLRKIQIKKLCASGIKE SPCVVYEDMSYSNRKTPCIPNQYQ
Source	in vitro E.coli expression system
Target Names	Cd7
Protein Names	Recommended name: T-cell antigen CD7 Alternative name(s): CD_antigen= CD7
Expression Region	24-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.