



Recombinant Rat T-cell surface glycoprotein CD8 beta chain (Cd8b)

Product Code	CSB-CF004967RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P05541
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	LLQTPSSLLVQTNQTAKMSCEAKTFPKGTTIYWLRELQDSNKNKHFEFLASRT STKGIKYGERVKKNMTLNLSFNSTLPFLKIMDVKPEDSGFYFCAMVGGSPMVVFGT GTKLTVVDVLPPTTAPTKKTTLLKKKQCPTPHPKTQKGLTCGLITLSLLVACILVLLV SLSVAIHFHCMRRRRARIHFMKQFHK
Source	in vitro E.coli expression system
Target Names	Cd8b
Protein Names	Recommended name: T-cell surface glycoprotein CD8 beta chain Alternative name(s): CD8 antigen 37 kDa chain OX-8 membrane antigen CD_antigen=CD8b
Expression Region	22-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.