



Recombinant Mouse Cytotoxic T-lymphocyte protein 4 (Ctla4)

Product Code	CSB-CF006163MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P09793
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	EAIQVTQPSVVLASSHGVASFPCEYSPSHNTDEV RVTVLRQTNDQMTEVCATT FTEKNTVGF LDY PFCSGTFN ES RVNLT IQGLRAVDTGLYLCKVELMYPPPYFV GMGNGTQIYVIDPEPCPDSDFLLWILVAVSLGLFFYSFLVSAVSLSKMLKKRSP LTTGVYVKMPPEPECEKQFQPYFIPIN
Source	in vitro E.coli expression system
Target Names	Ctla4
Protein Names	Recommended name: Cytotoxic T-lymphocyte protein 4 Alternative name(s): Cytotoxic T-lymphocyte-associated antigen 4 Short name= CTLA-4 CD_antigen= CD152
Expression Region	36-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.