



Recombinant Human Myelin protein zero-like protein 2 (MPZL2)

Product Code	CSB-CF014776HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O60487
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	VEIYTSRVLEAVNGTDARLKCTFSSFAPVGDALVTWNFRPLDGGPEQFVFYY HIDPFQPMGRFKDRVSWDGNPERYDASILLWKLQFDDNGTYTCQVKNPPDV DGVIGEIRLSVVHTVRFSEIHFLALAIGSACALMIIIVVVVLFQHYRKKRWAERA HKVVEIKSKEEERLNQEKKVSVYLEDTD
Source	in vitro E.coli expression system
Target Names	MPZL2
Protein Names	Recommended name: Myelin protein zero-like protein 2 Alternative name(s): Epithelial V-like antigen 1
Expression Region	27-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.