



Recombinant Mouse B-cell antigen receptor complex-associated protein alpha chain (Cd79a)

Product Code	CSB-CF004957MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11911
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	LRVEGGPPSLTVNLGEEARLTCENNGRNPNTWWFSLQSNITWPPVPLGPGQ GTTGQLFFPEVNKNHRGLYWCQVIENNILKRSCGTYLRVRNPVPRPFLDMGE GTKNRIITAEGIILLFCVVPGTLLLFRKRWQNEKFGVDMPPDDYEDENLYEGLN LDDCSMYEDISRGLQGTYQDVGNLHIGDAQLEKP
Source	in vitro E.coli expression system
Target Names	Cd79a
Protein Names	Recommended name: B-cell antigen receptor complex-associated protein alpha chain Alternative name(s): Ig-alpha MB-1 membrane glycoprotein Membrane-bound immunoglobulin-associated protein Surface IgM-associated protein CD_antigen= CD79a
Expression Region	29-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.