



Recombinant Mouse CD70 antigen (Cd70)

Product Code	CSB-CF004954MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O55237
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MPEEGRPCPWVRWSGTAFQRQWPWLLLVVFITVFCCWFHCSGLLSKQQQRL LEHPEPHTAELQLNLTVPKDPTRLRWGAGPALGRSFTHGPELEEGHLRIHQD GLYRLHIQVTLANCSSPGSTLQHRATLAVGICSPAAGHISLLRGRFGQDCTVAL QRLTYLVHGDVLCNTLTLPLPSRNADETFFGVQWICP
Source	in vitro E.coli expression system
Target Names	Cd70
Protein Names	Recommended name: CD70 antigen Alternative name(s): CD27 ligand Short name= CD27-L Tumor necrosis factor ligand superfamily member 7 CD_antigen= CD70
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.