



Recombinant Rabbit T-cell-specific surface glycoprotein CD28 (CD28)

Product Code	CSB-CF004913RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P42069
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	NKILVKQSPMLVVNNNEVNLSCKYTYNLFKSKEFRASLYKGADSAVEVCCVNGN FSHPHQFHSTTGFNCDGKLGNETVTFYLKNLYVNQTDIYFCKIEVMYPPPYLD NEKSNGTIIHVKEQHFCAHPSPKSSTLFWVLVVVGAVLAFYSMLVTVALFSC WMKSKKNRLLQSDYMNMTPRRPGPTRKHYPYAPARDFAAAYRS
Source	in vitro E.coli expression system
Target Names	CD28
Protein Names	Recommended name: T-cell-specific surface glycoprotein CD28 Alternative name(s): CD_antigen= CD28
Expression Region	20-221
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.