



Recombinant Mouse B-cell antigen receptor complex-associated protein beta chain (Cd79b)

Product Code	CSB-CF004958MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P15530
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	VPAMTSSDLPLNFGGSPCSQIWQHPRFAAKKRSSMVKFHCYTNHSGALTWF RKRGSQQPQELVSEEGRIVQTQNGSVYTLTIQNIQYEDNGIYFCKQKCDSANH NVTDSCGTELLVLGFSTLDQLKRRNTLKDGIILIQTLIIILFIIVPIFLLLDKDDGKA GMEEDHTYEGLNIDQTATYEDIVTLRTGEVKWSVGEHPGQE
Source	in vitro E.coli expression system
Target Names	Cd79b
Protein Names	Recommended name: B-cell antigen receptor complex-associated protein beta chain Alternative name(s): B-cell-specific glycoprotein B29 Ig-beta Immunoglobulin-associated B29 protein CD_antigen= CD79b
Expression Region	26-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.