



Recombinant Danio rerio Apoptosis-related protein 3 (apr3)

Product Code	CSB-CF003836DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3KNS9
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	QDTDAQLCQMCEGTIRHDSFVWSFCITKGYVKGHCCKFNNTSDVDTIIGLDLS NCSISHVEHLYNSSTALIIDLNNPISNLSQDYVFQGFSQLTQLLLPSKLECPGGR ASWEKVEVKSITRICEGQKNACNQTVQMPLVCPENSLCSPYGPFFECCLN NFHGYKCMRQGEFPLVKVLGILTASTVVVSSVLWFTQRRKVKNT
Source	in vitro E.coli expression system
Target Names	atraid
Protein Names	Recommended name: Apoptosis-related protein 3 Short name= APR-3
Expression Region	27-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.