



Recombinant Human Fms-related tyrosine kinase 3 ligand (FLT3LG)

Product Code	CSB-CF008734HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49771
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	TQDCSFQHSPISSDFAVKIRELSDYLLQDYPVTVASNLQDEELCGGLWRLVLA QRWMERLKTAVAGSKMQGLLERVNTTEIHFVTKCAFQPPPSCLRFVQTNISRLLQ ETSEQLVALKPWITRQNFRCLELQCQPDSSTLPPPWSRPLEATAPTAPQPP LLLLLLLPVGLLLLAAAWCLHWQRTRRRTPRPGEQVPPVPSPQDLLLVEH
Source	in vitro E.coli expression system
Target Names	FLT3LG
Protein Names	Recommended name: Fms-related tyrosine kinase 3 ligand Short name= Flt3 ligand Short name= Flt3L Alternative name(s): SL cytokine
Expression Region	27-235
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.