



Recombinant Human Golgi SNAP receptor complex member 2 (GOSR2)

Product Code	CSB-CF009678HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14653
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MDPLFQQTHKQVHEIQSCMGRLETADKQSVHIVENEIQASIDQIFSRRLERLEILS SKEPPNKRQNARLRVDQLKYDVQHLQTALRNFAQHRRHAREQQERQREELLS RTFTTNDSDTTIPMDES LQFNSSLQKVHNGMDDLILDGHNILDGLRTQRLTLKG TQKKILDIANMLGLSNTVMRLIEKRAFQDKYFMIGGMLLTCVVMFLVVQYLT
Source	in vitro E.coli expression system
Target Names	GOSR2
Protein Names	Recommended name: Golgi SNAP receptor complex member 2 Alternative name(s): 27 kDa Golgi SNARE protein Membrin
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.