



Recombinant Human Receptor-binding cancer antigen expressed on SiSo cells (EBAG9)

Product Code	CSB-CF007355HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O00559
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAITQFRLFKFCTCLATVFSFLKRLICRSGRGRKLSGDQITLPTTVDYSSVPKQT DVEEWTSWDEDAPTSVKIEGGNGNVATQQNSLEQLEPDYFKDMTPTIRKTQKI VIKKREPLNFGIPDGSTGFSSRLAATQDLPFQHSSELGDLDTWQENTNAWEE EEDAAWQAEVLRQQKLADREKRAAEQQRKKMEKEAQRLMKKEQNKIGVKL S
Source	in vitro E.coli expression system
Target Names	EBAG9
Protein Names	Recommended name: Receptor-binding cancer antigen expressed on SiSo cells Alternative name(s): Cancer-associated surface antigen RCAS1 Estrogen receptor-binding fragment-associated gene 9 protein
Expression Region	1-213
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.