



Recombinant Human Immunoglobulin superfamily member 6 (IGSF6)

Product Code	CSB-CF011564HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95976
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	CTLSVTQPWYLEVDYTHEAVTIKCTFSATGCPSEQPTCLWFRYGAHQPENLCL DGCKSEADKFTVREALKENQVSLTVNRVTSNDSAIYICGIAFPSVPEARAKQTG GGTTLVVREIKLLSKELRSFLTALVSLLSVYVTGVCVAFILLSKSKSNPLRNKEIK EDSQKKKSARRIFQEIAQELYHKRHVETNQQSEKDNNTYENRRVLSNYERP
Source	in vitro E.coli expression system
Target Names	IGSF6
Protein Names	Recommended name: Immunoglobulin superfamily member 6 Short name= IgSF6 Alternative name(s): Protein DORA
Expression Region	28-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.