



Recombinant Human Vesicle-trafficking protein SEC22b (SEC22B)

Product Code	CSB-CF020945HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O75396
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	VLLTMIARVADGLPLAASMQEDEQSGRDLQQYQSQAKQLFRKLNEQSPTRCT LEAGAMTFHYIIEQGVCDLVLCEAAFPKTLAFAYLEDLHSEFDEQHGKKVPTVS RPYSFIEFDTFIQKTKKLYIDSCARRNLGSINTELQDVQRIMVANIEEVLQRGEA LSALDSKANNLSSLSKKYRQDAKYLMHSTYAKLAAVAVFFIMLIVYVRFWWL
Source	in vitro E.coli expression system
Target Names	SEC22B
Protein Names	Recommended name: Vesicle-trafficking protein SEC22b Alternative name(s): ER-Golgi SNARE of 24 kDa Short name= ERS-24 Short name= ERS24 SEC22 vesicle-trafficking protein homolog B SEC22 vesicle-trafficking protein-like 1
Expression Region	2-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.