



Recombinant Rat NKG2-D type II integral membrane protein (Klrk1)

Product Code	CSB-CF012474RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O70215
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MSKCHNYDLKPAKWDTSQEHQKQRSALPTSRPGengiIRRRSSIEELKISPLFV VRVLVAAMTIRFTVITLTWLAVFITLLCNKEVSVSSREGYCGPCPNDWICHRNN CYQFFNENKAWNQSQASCLSQNSSLLKIYSKEEQDFLKLVKSYHWMGLVQSP ANGSWQWEDGSSLSPNELTLVKTPSGTCAVYGSSFKAYTEDCSNPNTYICMK RAV
Source	in vitro E.coli expression system
Target Names	Klrk1
Protein Names	Recommended name: NKG2-D type II integral membrane protein Alternative name(s): Killer cell lectin-like receptor subfamily K member 1 NK cell receptor D NK lectin-like receptor Short name= NKLLR NKG2-D-activating NK receptor NKR-P2 CD_ant
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.