



Recombinant Macaca fascicularis NKG2-D type II integral membrane protein (KLRK1)

Product Code	CSB-CF012474MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61252
Product Type	Transmembrane Protein
Immunogen Species	Macaca fascicularis (Crab-eating macaque) (Cynomolgus monkey)
Sequence	MGWIRGRRPRHNLEMSEFHNYKLGLAKSDFSTRCQKQRCPVIKSKCRENASP LFFCCFIAMGIRFIIMVTIWSAVFLNSLFNQEVQIPLTESYCGPCPKNWICYK NNCYQFFNESKNWYESQASCMSQNASLLKVYSKEDQDLLKLVKSYHWMGLV HIPTNGSWQWEDGSILSPNLLTIEMQKGDALYASSFKGYIENC SIPNTYICMQ RTV
Source	in vitro E.coli expression system
Target Names	KLRK1
Protein Names	Recommended name: NKG2-D type II integral membrane protein Alternative name(s): Killer cell lectin-like receptor subfamily K member 1 NK cell receptor D NKG2-D-activating NK receptor CD_antigen= CD314
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.