



Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 3 (CEACAM3)

Product Code	CSB-CF005163HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P40198
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	KLTIESMPLSVAEGKEVLLLVHNLQPQLFGYSWYKGERVDGNSLIVGYVIGTQQ ATPGAAYSGRETIYTNASLLIQNVTDNDIGFYTLQVIKSDLVNEEATGQFHVYQ ENAPGLPVGAVAGIVTGVLVGVVALVAALVCFLLAKTGRTSIQRDLKEQQPQAL APGRGPSHSSAFMSPLSTAQAPLPNPRTAASIYEELLKHDTNIYCRMDHKAE VAS
Source	in vitro E.coli expression system
Target Names	CEACAM3
Protein Names	Recommended name: Carcinoembryonic antigen-related cell adhesion molecule 3 Alternative name(s): Carcinoembryonic antigen CGM1 CD_antigen= CD66d
Expression Region	35-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.