



Recombinant Dog T-cell surface glycoprotein CD8 alpha chain (CD8A)

Product Code	CSB-CF004966DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P33706
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	SGPSRFRMTPPKVVGQLHAQVELQCQVLLSTAAPGCSWLYQRNEPAARPVFL MYISQSRAKPAEGLDTKHISGQKKTSTYSLTLRFRKEDEGYFCSVLSNSIL YFSPFVPVFLPVKPPTTPAPRPPTTRAPTNAKPVSPRGETCRPAAGSAVKTSGL LDFACEIYWAPLAGTCAVLLLSLVITIICNHRNRRRVCKCPRPVVRPGGKPSPS EKYV
Source	in vitro E.coli expression system
Target Names	CD8A
Protein Names	Recommended name: T-cell surface glycoprotein CD8 alpha chain Alternative name(s): CD_antigen= CD8a
Expression Region	22-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.