



Recombinant Cat T-cell surface glycoprotein CD8 alpha chain (CD8A)

Product Code	CSB-CF004966CA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41688
Product Type	Transmembrane Protein
Immunogen Species	Felis catus (Cat) (Felis silvestris catus)
Sequence	AGPSPFRLSPVRVEGRLGQRVELQCEVLLSSAAPGCTWLFQKNEPAARPIFLA YLSRSRTKLAEEELDPKQISGQRIQDTLYSLTLHRFRKEEEGYFCSVVSNSVLY FSAFVPVFLPVKPTTTPAPRPPTQAPITTSQRVSLRPGTCQPSAGSTVEASGL DLSCDIYWAPLAGTCAFLLSLVITVICNHRNRRRVCKCPRPVVRAGGKPSPS ERYV
Source	in vitro E.coli expression system
Target Names	CD8A
Protein Names	Recommended name: T-cell surface glycoprotein CD8 alpha chain Alternative name(s): CD_antigen= CD8a
Expression Region	22-239
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.