



Recombinant Human Kunitz-type protease inhibitor 2 (SPINT2)

Product Code	CSB-CF022585HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O43291
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	ADRERSIHDFCLVSKVVGRCRASMPRWWYNVTDGSCQLFVYGGCDGNSNN YLTKEECLKKCATVTENATGDLATSRNAADSSVPSAPRRQDSEDHSSDMFNY EEYCTANAVTGPCRASFPRWYFDVERNSCNNFIYGGCRGNKNSYRSEEACM LRCFRQQENPPLPLGSKVVVLAGLFVMVLILFLGASMVYLIRVARRNQERALRT VWSSGDDKEQLVKNTYVL
Source	in vitro E.coli expression system
Target Names	SPINT2
Protein Names	Recommended name: Kunitz-type protease inhibitor 2 Alternative name(s): Hepatocyte growth factor activator inhibitor type 2 Short name= HAI-2 Placental bikunin
Expression Region	28-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.