



Recombinant Human NKG2-C type II integral membrane protein (KLRC2)

Product Code	CSB-CF012467HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P26717
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSKQRGTFSEVSLAQDPKRQQRKPKGNKSSISGTEQEIFQVELNLQNPSLNH QGIDKIYDCQGLLPPEKLTAEVLGIICIVLMATVLKTIVLIPFLEQNNSSPNTRTQ KARHCGHCPEEWITYSNSCYIGKERRTWEESLLACTSKNSSLLSIDNEEEMK FLASILPSSWIGVFRNSSHHPWVTINGLAFKHKIKDSDNAELNCAVLQVNRLLKS AQCGSSMIYHCKHKL
Source	in vitro E.coli expression system
Target Names	KLRC2
Protein Names	Recommended name: NKG2-C type II integral membrane protein Alternative name(s): CD159 antigen-like family member C NK cell receptor C NKG2-C-activating NK receptor CD_antigen= CD159c
Expression Region	1-231
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.