



Recombinant Human Myelin protein zero-like protein 1 (MPZL1)

Product Code	CSB-CF014775HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95297
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	SALEVYTPKEIFVANGTQGKLTCKFKSTSTTGGLTSVSWSFQPEGADTTVSFF HYSQGQVYLGNYPPFKDRISWAGDLDDKKDASINIENMQFIHNGTYICDVKNPP DIVVQPGHIRLYVVEKENLPVFPVWVVVGIVTAVVLGLTLLISMILAVLYRRKNS KRDTGCGSTSESLSPVKQAPRKSPDTEGLVKSLPSGSHQGPVIYAQLDHSG GHHSDKINKSESVVYADIRKN
Source	in vitro E.coli expression system
Target Names	MPZL1
Protein Names	Recommended name: Myelin protein zero-like protein 1 Alternative name(s): Protein zero-related
Expression Region	36-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.