



Recombinant Rat C-type lectin domain family 9 member A (Clec9a)

Product Code	CSB-CF005542RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	D4AD02
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MHEEEIYTSLQWDIPTSEASQKCPSSLKCPGTWCIVTVISCVVVCVGLLAASIFLG IKFSQVSSLVMEQRERLIRQDTALLNLTEWQRNHTLQLKSCQASLQRSLRSGS NCNPPNWIQNGKSCYYAFDRWETWNNSSKSKLKEGDSLLQIDSKEEMEFI NLSIWKLKGGYEWVGVFQDGPSGSWFEDGSSPLSDLLPTDRQLSASQIC GYLKDHTLISDNCSNWKYFICEKKAFGSCI
Source	in vitro E.coli expression system
Target Names	Clec9a
Protein Names	Recommended name: C-type lectin domain family 9 member A
Expression Region	1-241
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.