



Recombinant Dog Interleukin-2 receptor subunit alpha (IL2RA)

Product Code	CSB-CF011649DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O62802
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	DLCDDDDPPNLKHATFKALTYKTGTVLNCDCERGFRRISSYMHCTGNSSHASW ENKCRCKSVSPENRKGKVTTKPEEQKGENPTEMQSQTTPMDEVLDLVGHCRE PPPWEHENSRIYHFVVGQTLHYQCMQGFTALHRGPAKSICKTIFGKTRWTQ PPLKCISESQFPDDEELQASTDAPAGRDTSSPFITTSTPDFHKHTEVATTMESF IFTTEYQIAVASCVLLISIVLLSGLTWQRRRRRSRTI
Source	in vitro E.coli expression system
Target Names	IL2RA
Protein Names	Recommended name: Interleukin-2 receptor subunit alpha Short name= IL-2 receptor subunit alpha Short name= IL-2-RA Short name= IL-2R subunit alpha Short name= IL2-RA Alternative name(s): TAC antigen p55 CD_antigen= CD25
Expression Region	22-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.