



Recombinant Rat OX-2 membrane glycoprotein (Cd200)

Product Code	CSB-CF004895RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04218
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	QVEVVTQDERKLLHTTASLRCSLKTQEPLIVTWQKKKAVGPENMVTYSKAHG VVIQPTYKDRINITELGLLNTSITFWNTTLDDEGCYMCLFNMFGSGKVSGTACL TLYVQPIVHLHYNYFEDHLNITCSATARPAPAIWKG TGSGIENSTESHSHSNG TTSVTSILRVKDPKTQVGKEVICQVLYLGNVIDYKQSLDKGFWSVPLLLSIVSL VILLVLISILLYWKRHRNQERGESQGMQRMK
Source	in vitro E.coli expression system
Target Names	Cd200
Protein Names	Recommended name: OX-2 membrane glycoprotein Alternative name(s): MRC OX-2 antigen CD_antigen= CD200
Expression Region	31-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.