



Recombinant Mouse Fc receptor-like protein 6 (Fcrl6)

Product Code	CSB-CF008559MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A1YIY0
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	QELFPNPELTEFTNSETMDVILKCTIKVDPKNPTLQLFYTFYKNNHVIQDRSPHS VFSAEAKEENSGLYQCMVDTEGLIQKKSGYLDIQFWTPVSHPVLTQLQHEATN LAVGDKVEFLCEAHQGSPLPIFYSFYINGEILGKPLAPSGRAASLLASVKAEWST KNYSCEAKNNISREISELKKFPLVVSGTAWIKSNMLPIWLPASLLGGMVIAAVVL MYFFKPKCKKHARPETPTLKEPDSFLYVSVDNQRYSK
Source	in vitro E.coli expression system
Target Names	Fcrl6
Protein Names	Recommended name: Fc receptor-like protein 6 Short name= FcR-like protein 6 Short name= FcRL6
Expression Region	17-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.