



Recombinant Mouse Tumor necrosis factor receptor superfamily member 4 (Tnfrsf4)

Product Code	CSB-CF023981MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P47741
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	VTARRLNCVKHTYPSGHKCCRECPGHGMVSRCDHTRDTLCHPCETGFYNE AVNYDTCKQCTQCNHRSGSELKQNCTPTQDTVCRCPGTQPRQDSGYKLGV DCVPCPPGHFSPGNNQACKPWTNCTLSGKQTRHPASDSLDAVCEDRSLLATL LWETQRPTFRPTTVQSTTVWPRTSELPSPPTLVTEGPAFAVLLGLGLLLAP LTVLLALYLLRKAWRLPNTPKPCWGNFSFRTPIQEEHTDAHFTLAKI
Source	in vitro E.coli expression system
Target Names	Tnfrsf4
Protein Names	Recommended name: Tumor necrosis factor receptor superfamily member 4 Alternative name(s): OX40 antigen OX40L receptor CD_antigen= CD134
Expression Region	20-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.