



Recombinant Cat Interleukin-2 receptor subunit alpha (IL2RA)

Product Code	CSB-CF011649CA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41690
Product Type	Transmembrane Protein
Immunogen Species	Felis catus (Cat) (Felis silvestris catus)
Sequence	ELCDENPPDIQHATFKALTYKTGTM LNCECKKGFRRISNGSAFMLCAGNSSHS SWENQRCRCISTSPRATDGQIIPKPEEQKGKSPMGMQSQMQPTDQVNLPGHC REPPPWEHEDSRRIYHFVVGQTVHYQCMQGFRA LKRGPAKSVCKTTTCGKAT WTQ PRLQCLSERSDGWFPDDEEPQASTDAALGSDTSCPSITASTTDFQRHTE VAMTTESFVFTTEYQIAVAGCVLLLISIVLLSGLTWQRRWRKSRRTI
Source	in vitro E.coli expression system
Target Names	IL2RA
Protein Names	Recommended name: Interleukin-2 receptor subunit alpha Short name= IL-2 receptor subunit alpha Short name= IL-2-RA Short name= IL-2R subunit alpha Short name= IL2-RA Alternative name(s): CD_antigen= CD25
Expression Region	22-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.