



Recombinant *Saccharomyces cerevisiae* Syntaxin-8 (SYN8)

Product Code	CSB-CF022904SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P31377
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MDVLKLG YELDQLSDLVEERTRLVSVLKLAPTSNDNVTLKRQLGSILELLQKCA PNDELISRYNTILDKIPDTAVDKELYRFQQQVARNTDEVSKESLKKVRFKNDDE LTVMYKDDDEQDEESPLPSTHTPYKDEPLQSQLQSQSQPPQPMVSNQEL FINQQQLLEQDShLGALSQSIGRTHDISLDLNNEIVSQNDSLLVDLENLIDNNG RNLNRASRSMHGFNNSRFKDNNGNCVILVLIVVLLLLLLVL
Source	in vitro E.coli expression system
Target Names	SYN8
Protein Names	Recommended name: Syntaxin-8 Alternative name(s): SNARE protein related to mammalian syntaxin 8 ULP1-interacting protein 2
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.