



Recombinant Bovine Tissue factor (F3)

Product Code	CSB-CF007928BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30931
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	TDVVVAYNITWKSTNFKTILEWEPKPINHVYTVQISPR LGNWKNKCFYTTNTEC DVTDEIVKNVRETYLARVLSYPADTSSSTVEPPFTNSPEFTPYLETNLGQPTIQ SFEQVGTKLNVTVQDARTLVRANSAFLSLRDVFGKDLNYTLYYWKASSTGKKK ATTNTNGFLIDVDKGENYCFHVQAVILSRRVNQKSPESPIKCTSHEKVLSTELF FIIGTVMLVIII FIVVLSVSLHKCRKVRAERSGKENTPLNAA
Source	in vitro E.coli expression system
Target Names	F3
Protein Names	Recommended name: Tissue factor Short name= TF Alternative name(s): Coagulation factor III CD_antigen= CD142
Expression Region	36-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.