



Recombinant Danio rerio Malectin (mlec)

Product Code	CSB-CF014620DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A9C3P0
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	ESGAQSLAERVIWAVNAGGDTHTDVHGIQFKKDPLEGKVGKASDYGVRLPILRSSPEDQILYQTERYNEDTFGYEVPPIREEGDYILVMKYAEVYFAQSQQKVFVDVRLNGHVVKDLDFDRVGHSTAHDIVPFSIKRGKLSVHGEVSTFNGKLTVEFVKGYDNPICALYVMKGKLEDVPKLQPHPGLEKREEEEEEEEEEGEGPEGEKKSASTSPKNPVRSGPRTPNPYATDNSSLMFPILVAFGVFIPTLFCLCRL
Source	in vitro E.coli expression system
Target Names	mlec
Protein Names	Recommended name: Malectin
Expression Region	27-285
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.